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

Histogram 8

- Histogram 8 displays events from region A.
- Viable events are included within Region J. The 7-AAD Viability Dye positive events are excluded from Region J.

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
Histogram 1

- Histogram 1 displays viable events from region J (see Histogram 8). If not using 7-AAD, Histogram 1 displays all events minus all Stem-Count Fluorospheres (i.e. « -H » see Histogram 5).

NOTE:

- Position Region A to include all CD45⁺ events (leukocytes) while excluding CD45⁻ events. Position Region E to include only lymphocytes (bright CD45, low Side Scatter).
- Region A is intended to serve as an appropriate denominator for viable WBC in the calculation of the percentage of CD34⁺ HPC. Once the analysis is done, the absolute count of viable leukocytes is given on the statistic printout related to Region A.

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
Histogram 2


- Histogram 2 displays events from region A and J. If not using 7-AAD, Histogram 2 displays events from region A.
- Adjust Region B to surround CD34⁺ HPC, including CD34^{dim} events. The upper limit along the y-axis may be set to include CD34⁺ events with low to intermediate Side Scatter for CD34⁺ HPC.

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
Histogram 3

- Histogram 3 displays events from A and B and J. If not using 7-AAD, Histogram 2 displays events from region A and B.
- Adjust Region C to include cells forming a cluster with characteristic $CD34^+$ HPC (i.e. low Side Scatter and low to intermediate CD45 staining). Note the characteristic shape of Region C. It is the cluster of events that determines where Region C is centered. Brightly FITC stained events (platelet aggregates and mature monomyeloid cells) must be excluded in the setting of this region.

Histogram 4

- Histogram 4 displays events from A and B and C and J. If not using 7-AAD, Histogram 2 displays events from region A and B and C.
- Lymph / Blast Region D identifies a cluster of events meeting all the fluorescence and light scatter criteria of ISHAGE Guidelines for $CD34^+$ HPC. Note the characteristic shape of Region D. If present, platelet aggregates that stain weakly must be excluded to the left of Region D.
- Once the analysis is complete, the absolute count of viable $CD34^+$ HPC is given on the statistic printout related to the $CD34^+$ Stem-Trol Control Cells or $CD34^+$ HPC (Region D).

Histogram 5 displays all events.


NOTE:

This histogram is useful to visualize the lower limit of CD45 expression within the $CD34^+$ events.

- Set Quadstat Region I2 to enclose all $CD45^+CD34^+$ events. Set the left boundary of Quadstat Region I2 to include all $CD34^+CD45^{dim}$ events. Verify on Histogram 1 that the lower limit of Region A includes all $CD34^+$ events (use Histogram 5 as a guide).
- **Amorphous** Region H is drawn to include all Stem-Count Fluorospheres. Region H must include the last channel of the FL1 Log (on the right) and FL2 Log scales (on the top).
- Set Region K (named "listgate -") to exclude most of the double negative events from the acquisition.

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

Histogram 6

- Histogram 6 displays events from E and J regions. If not using 7-AAD, Histogram 2 displays events from region E.
- Verify whether the Forward and Side Scatter gain parameters of the flow cytometer are optimally set for the processed sample. If necessary, adjust the Forward Scatter voltage / gain so that the smallest lymphocytes scatter in the middle-left part of the histogram. Adjust the Forward Scatter discriminator to ensure that even the smallest lymphocytes scatter above the discriminator.

Histogram 7

- Histogram 7 displays events from region H.
- Region G encloses the Stem-Count Fluorospheres singlet population. Check that the fluorosphere singlets accumulate homogeneously and constantly over time. Region G can be labeled as $\leq \text{CAL} \leq$ to allow automatic calculation of absolute numbers of CD34^+ HPC. Type the correct assayed concentration of the current lot of Stem-Count Fluorospheres (refer to the next section of this insert or to the instrument manual for further details).

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